

Venom-D Prototype


By Hemanth Sabbella

Introduction

- The idea to detect various snake species venom in human blood to do a one-to-one mapping of venom and the snake species using Venom Detection Kit
- Currently, we are detecting a specific venom in gold nanoparticles sample using a low-cost handheld device.
- Keynote: Gold nanoparticles has the an absorbance peak at 520 nm wavelength and when added with a specific venom there will be a wavelength peak shift.

Beer-Lambert's Law

"The Beer-Lambert law states that the quantity of light absorbed by a substance dissolved in a fully transmitting solvent is directly proportional to the concentration of the substance and the path length of the light through the solution."



Transmittance:

$$T = I / I_0$$

$$\% \text{ Transmittance: } \%T = 100 T$$

Absorbance:


$$A = \log_{10} (I_0/I)$$

$$A = \log_{10} (1/T) = -\log_{10} (T)$$


$$A = \log_{10} (100/\%T)$$

$$A = 2 - \log_{10} (\%T)$$

Conventional Spectrophotometer




What's in the box?




- SMD LED as transmitting sources.
- Photodetector array IC board
- Arduino Nano microcontroller.

Transmitter Board

- Using narrow-band, dominant wavelength LEDs as light sources for a specific wavelength absorbance.
- 6-point light sources.






Spectral Response: Wavelength vs Luminous Intensity




Receiver Board


- TCS3200 - Light(irradiance) to frequency converter.
- 8 x 8 array of photodiodes.
 - 16 blue filters, 16 green filters, 16 red filters and 16 clear filters.
- The output is a square wave (50% duty cycle) with frequency directly proportional to luminous intensity (irradiance).




Cheat Sheet: 3D modeling and printing





Modeling



Printing: Calibrating the printer



Slicing



Printer Configuration


Print Setup




Note: Need to use generate support option for some models; filament won't stick in air.

Signal conditioning & Calibration

- Zero referencing.
 - Luminous intensity mismatch.
 - For every measurement, a initial referencing to start with.



- Range Mapping
 - Same scale: Relative change in one reading has to be mapped to change in other.



THANK YOU!